



WAC Mouse mAb

Catalog No	YP-mAb-18599
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human WAC (NP_057712.2).
Specificity	
Formulation	
Source	
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	11167P7Rik;BA48B24;BA48B24.1;BM 16;FLJ3129;KIAA1844;MGC1753;OTTHUMP1938;OTTHUMP19382;PRO1741;wac;WAC_HUMAN;WW domain containing adapter protein with coiled coil;WW domain containing adaptor with coiled coil;WW domain-containing adapter protein with coiled-coil;Wwp4;DESSH;BM-016
Observed Band	71kDa
Calculated Molecular Weight	71kD(Calculated).
Cell Pathway	
Tissue Specificity	
Function	
Background	The protein encoded by this gene contains a WW domain, which is a protein module found in a wide range of signaling proteins. This domain mediates protein-protein interactions and binds proteins containing short linear peptide motifs that are proline-rich or contain at least one proline. This gene product shares 94% sequence identity with the WAC protein in mouse, however, its exact function is not known. Alternative splicing results in multiple transcript



variants


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of various lysates using WAC Mouse pAb (A17703) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 60s