



SAPS3 Mouse mAb

Catalog No	YP-mAb-18606
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PPP6R3 C11orf23 KIAA1558 PP6R3 SAPL SAPS3
Immunogen	Synthesized peptide derived from human SAPS3
Specificity	This antibody detects endogenous levels of SAPS3 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SAPLa;SAPL;PPP6R3;PP6R3;KIAA1558;PPP6R3 C11orf23 KIAA1558 PP6R3 SAPL SAPS3;Serine/threonine-protein phosphatase 6 regulatory subunit 3 (SAPS domain family member 3) (Sporulation-induced transcript 4-associated protein SAPL)
Observed Band	130-140 kDa
Calculated Molecular Weight	873 aa, 98 kDa
Cell Pathway	Cytoplasm. Nucleus.
Tissue Specificity	Expressed in skeletal muscle, placenta, heart, pancreas, testis, brain, lung, liver, kidney, spleen, thymus, prostate, small intestine, colon and leukocytes.
Function	Regulatory subunit of protein phosphatase 6 (PP6). May function as a scaffolding PP6 subunit. May have an important role in maintaining immune self-tolerance.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
