



EDIL3 Mouse mAb

Catalog No	YP-mAb-18644
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EDIL3 DEL1
Immunogen	Synthesized peptide derived from human EDIL3
Specificity	This antibody detects endogenous levels of EDIL3 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DEL1;Developmental endothelial locus 1;Developmentally regulated endothelial cell locus 1 protein;Developmentally-regulated endothelial cell locus 1 protein; Edil3;EDIL3_HUMAN;EGF like repeat and discoidin I like domain containing protein 3;EGF like repeats and discoidin I like domains 3;EGF like repeats and discoidin I like domains protein 3;EGF-like repeat and discoidin I-like domain-containing protein 3;Integrin binding protein DEL1;Integrin-binding protein DEL1; EDIL3 DEL1;EGF-like repeat and discoidin I-like domain-containing protein 3 (Developmentally-regulated endothelial cell locus 1 protein) (Integrin-binding protein DEL1)
Observed Band	53kD
Calculated Molecular Weight	52-55 kDa; 54kD(Calculated).
Cell Pathway	Secreted.
Tissue Specificity	
Function	Promotes adhesion of endothelial cells through interaction with the alpha-v/beta-3 integrin receptor. Inhibits formation of vascular-like structures. May be involved in regulation of vascular morphogenesis of remodeling in embryonic development.



Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images