



RPUSD2 Mouse mAb

Catalog No	YP-mAb-18647
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPUSD2 C15orf19
Immunogen	Synthesized peptide derived from human RPUSD2
Specificity	This antibody detects endogenous levels of RPUSD2 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RNA pseudouridylate synthase domain-containing protein 2;Pseudouridylate synthase RPUSD2;EC:5.4.99.-;C18B11;C15orf19;RPUSD2 C15orf19
Observed Band	80-90 kDa
Calculated Molecular Weight	545 aa, 61 kDa
Cell Pathway	
Tissue Specificity	
Function	Pseudouridine synthase that catalyzes pseudouridylation of mRNAs.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
