



Catalog No	YP-mAb-18648
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZCRB1
Immunogen	Synthesized peptide derived from human ZCRB1
Specificity	This antibody detects endogenous levels of ZCRB1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Zinc finger CCHC-type and RNA-binding motif-containing protein 1;SNRNP31; ZCRB1;Zinc finger CCHC-type and RNA-binding motif-containing protein 1 (U11/U12 small nuclear ribonucleoprotein 31 kDa protein) (U11/U12 snRNP 31 kDa protein) (U11/U12-31K)
Observed Band	25-30 kDa
Calculated Molecular Weight	217 aa, 25 kDa
Cell Pathway	Nucleus, nucleoplasm .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
