




HSP20/HSPB6 Mouse mAb

Catalog No	YP-mAb-18651
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EXOSC8 OIP2 RRP43
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 70-150 of human HSP20/HSPB6 (O14558)
Specificity	
Formulation	
Source	
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HEL55;Hsp20;PPP1R91;HSP20/HSPB6
Observed Band	17kDa
Calculated Molecular Weight	17kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	This locus encodes a heat shock protein. The encoded protein likely plays a role in smooth muscle relaxation.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of various lysates using HSP20/HSPB6 Mouse mAb (A9091) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 1s.