



VDAC2 Mouse mAb

Catalog No	YP-mAb-18670
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VDAC2
Immunogen	Synthesized peptide derived from human VDAC2
Specificity	This antibody detects endogenous levels of VDAC2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	POR;VDAC2;Voltage-dependent anion-selective channel protein 2 (VDAC-2) (hVDAC2) (Outer mitochondrial membrane protein porin 2)
Observed Band	30kDa
Calculated Molecular Weight	32kDa
Cell Pathway	Mitochondrion outer membrane . Membrane . May localize to non-mitochondrial membranes. .
Tissue Specificity	Expressed in erythrocytes (at protein level) (PubMed:27641616). Expressed in all tissues examined (PubMed:8420959).
Function	Forms a channel through the mitochondrial outer membrane that allows diffusion of small hydrophilic molecules (By similarity). The channel adopts an open conformation at low or zero membrane potential and a closed conformation at potentials above 30-40 mV (By similarity). The open state has a weak anion selectivity whereas the closed state is cation-selective (By similarity). Binds various lipids, including the sphingolipid ceramide, the phospholipid phosphatidylcholine, and the sterol cholesterol . Binding of ceramide promotes the mitochondrial outer membrane permeabilization (MOMP) apoptotic pathway .
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images