




CHCHD7 Mouse mAb

Catalog No	YP-mAb-18691
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-97 of human CHCHD7 (NP_077276.2)
Specificity	
Formulation	
Source	
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	COX23;CHCHD7
Observed Band	10kDa
Calculated Molecular Weight	10kDa
Cell Pathway	
Tissue Specificity	
Function	Acts as a transcriptional repressor .
Background	Predicted to be located in mitochondrial intermembrane space
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Mouse kidney, using CHCHD7 Mouse pAb (A18489) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 10s