




FBXO33 Mouse mAb

Catalog No	YP-mAb-18697
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 346-555 of human FBXO33 (NP_976046.1).
Specificity	
Formulation	
Source	
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F box only protein 33;F-box only protein 33;Fbx33;FBX33_HUMAN;FBXO33;BMND12;c14_5247
Observed Band	63kDa
Calculated Molecular Weight	7kD; 63kD(Calculated).
Cell Pathway	
Tissue Specificity	
Function	
Background	This locus represents an member of the F-box gene family. The encoded protein contains an F-box motif and a domain that might form a structure similar to a leucine-rich repeat found in placental RNase inhibitor. This locus may be associated with copy number variation of UGT2B17 (GeneID 7367), which has been associated with susceptibility to osteoporosis
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of various lysates using FBXO33 Mouse pAb (A14303) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 20s.