



ARL6IP4 Mouse mAb

Catalog No	YP-mAb-18698
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 313-402 of human ARL6IP4 (NP_057722.3).
Specificity	
Formulation	
Source	
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SR-25;SRp25;SFRS20;SRp37;ARL6IP4
Observed Band	25kDa
Calculated Molecular Weight	26kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Enables identical protein binding activity. Predicted to be involved in RNA splicing and mRNA processing. Predicted to be located in nucleus
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
