



PBX2 Mouse mAb

Catalog No	YP-mAb-18746
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PBX2 G17
Immunogen	Synthesized peptide derived from human PBX2
Specificity	This antibody detects endogenous levels of PBX2 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	G17;Homeobox protein PBX2;HOX12;PBX2;PBX2MHC;pre B cell leukemia homeobox 2;Protein G17;PBX2 G17;Pre-B-cell leukemia transcription factor 2 (Homeobox protein PBX2) (Protein G17)
Observed Band	50-60 kDa
Calculated Molecular Weight	430 aa, 46 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitously expressed.
Function	Transcriptional activator that binds the sequence 5'-ATCAATCAA-3'. Activates transcription of PF4 in complex with MEIS1.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
