



Calsequestrin 2 Mouse mAb

Catalog No	YP-mAb-18767
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CASQ2
Immunogen	Synthesized peptide derived from human Calsequestrin 2
Specificity	This antibody detects endogenous levels of Calsequestrin 2 at Human, Mouse, Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Calsequestrin 2;CASQ2;PDIB2;Calsequestrin-2 (Calsequestrin;cardiac muscle isoform); cardiac muscle isoform)
Observed Band	50 kDa
Calculated Molecular Weight	46 kDa
Cell Pathway	Sarcoplasmic reticulum lumen . This isoform of calsequestrin occurs in the sarcoplasmic reticulum's terminal cisternae luminal spaces of cardiac and slow skeletal muscle cells. .
Tissue Specificity	
Function	Calsequestrin is a high-capacity, moderate affinity, calcium-binding protein and thus acts as an internal calcium store in muscle. Calcium ions are bound by clusters of acidic residues at the protein surface, especially at the interface between subunits. Can bind around 60 Ca(2+) ions. Regulates the release of luminal Ca(2+) via the calcium release channel RYR2; this plays an important role in triggering muscle contraction. Plays a role in excitation-contraction coupling in the heart and in regulating the rate of heart beats.
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images