



C7orf58 Mouse mAb

Catalog No	YP-mAb-18770
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CPED1 C7orf58 UNQ9432/PRO34713
Immunogen	Synthesized peptide derived from human C7orf58
Specificity	This antibody detects endogenous levels of C7orf58 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cadherin-like and PC-esterase domain-containing protein 1;CPED1;UNQ9432/PRO34713;CPED1 C7orf58 UNQ9432/PRO34713
Observed Band	89 kDa
Calculated Molecular Weight	1026 aa, 117 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

