



ZNF395 Mouse mAb

Catalog No	YP-mAb-18776
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF395 HDBP2 PBF
Immunogen	Synthesized peptide derived from human ZNF395
Specificity	This antibody detects endogenous levels of ZNF395 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DKFZp434K121;HD gene regulatory region binding protein 2;HD gene regulatory region-binding protein 2;HD-regulating factor 2;HDBP-2;HDBP2;HDRF-2;Huntington disease gene regulatory region-binding protein 2;Huntington's disease gene regulatory region binding protein 2;Papillomavirus regulatory factor 1;papillomavirus regulatory factor PRF 1;Papillomavirus-binding factor;PBF;PRF-1;PRF1;PTTG1IP;Si-1-8-14;Zinc finger protein 395;ZN395_HUMAN;znf395;ZNF395 HDBP2 PBF;Zinc finger protein 395 (HD-regulating factor 2) (HDRF-2) (Huntington disease gene regulatory region-binding protein 2) (HD gene regulatory region-binding protein 2) (HDBP-2) (Papillomavirus regulatory factor 1) (PRF-1) (Papillomavirus-binding factor)
Observed Band	56kD
Calculated Molecular Weight	55kD(Calculated).
Cell Pathway	Cytoplasm. Nucleus. May shuttle between nucleus and cytoplasm.
Tissue Specificity	Widely expressed.
Function	Plays a role in papillomavirus genes transcription.
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images