



EAF2 Mouse mAb

Catalog No	YP-mAb-18782
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EAF2 TRAITS BM-040
Immunogen	Synthesized peptide derived from human EAF2
Specificity	This antibody detects endogenous levels of EAF2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BM 4;BM4;EAF 2;EAF2;EAF2_HUMAN;Ehrlich S II transcriptional activator factor;ELL associated factor 2;ELL-associated factor 2;Festa;Testosterone regulated apoptosis inducer and tumor suppressor;Testosterone regulated apoptosis inducer and tumor suppressor protein;Testosterone-regulated apoptosis inducer and tumor suppressor protein;TRAITS;U19;Uncharacterized bone marrow protein BM 4;Uncharacterized bone marrow protein BM4;EAF2 TRAITS BM-040;ELL-associated factor 2 (Testosterone-regulated apoptosis inducer and tumor suppressor protein)
Observed Band	29kD
Calculated Molecular Weight	29kD(Calculated).
Cell Pathway	Nucleus speckle .
Tissue Specificity	Expressed in heart, brain, placenta, lung, skeletal muscle, kidney, pancreas, spleen, prostate, testis, small intestine, colon, adrenal, bone marrow, lymph node, spinal gland, stomach, thyroid, trachea, thymus, liver and leukocytes.
Function	Acts as a transcriptional transactivator of TCEA1 elongation activity (By similarity). Acts as a transcriptional transactivator of ELL and ELL2 elongation activities. Potent inducer of apoptosis in prostatic and non-prostatic cell lines. Inhibits prostate tumor growth in vivo.



Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images