



MOXD1 Mouse mAb

Catalog No	YP-mAb-18786
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MOXD1 MOX UNQ2493/PRO5780
Immunogen	Synthesized peptide derived from human MOXD1
Specificity	This antibody detects endogenous levels of MOXD1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DBH like monooxygenase protein 1;DBH-like monooxygenase protein 1; dJ248E1.1;DKFZp564G22;Monooxygenase X;monooxygenase;DBH like 1;MOX; Moxd1;MOXD1_HUMAN;PRO578;MOXD1 MOX UNQ2493/PRO5780;DBH-like monooxygenase protein 1 (Monooxygenase X)
Observed Band	67kD
Calculated Molecular Weight	7kD(Calculated).
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type I membrane protein .
Tissue Specificity	Highly expressed in lung, kidney, brain and spinal cord.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

