



TRIM15 Mouse mAb

Catalog No	YP-mAb-18792
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRIM15 RNF93 ZNF178 ZNFB7
Immunogen	Synthesized peptide derived from human TRIM15
Specificity	This antibody detects endogenous levels of TRIM15 at Human
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RING finger protein 93;RNF93;TRIM15;tripartite motif containing 15;Zinc finger protein 178;Zinc finger protein B7;ZNF178;ZNFB7;TRIM15 RNF93 ZNF178 ZNFB7;Tripartite motif-containing protein 15 (RING finger protein 93) (Zinc finger protein 178) (Zinc finger protein B7)
Observed Band	52 kDa
Calculated Molecular Weight	465 aa, 52 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

