



MTHFD1 Mouse mAb

Catalog No	YP-mAb-18800
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MANEA
Immunogen	Synthesized peptide derived from human MANEA
Specificity	This antibody detects endogenous levels of MANEA at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CIMAH;MTHFC;MTHFD;MTHFD1;MANEA;Glycoprotein endo-alpha-1;2-mannosidase (Endo-alpha mannosidase) (Endomannosidase) (hEndo) (Mandaselin)
Observed Band	115kDa
Calculated Molecular Weight	102kDa
Cell Pathway	Highly expressed in the liver and kidney. Expressed at lower levels in muscle, pancreas, heart, placenta, lung and brain.
Tissue Specificity	
Function	Golgi apparatus membrane ; Single-pass type II membrane protein
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images