



CRIP2 Mouse mAb

Catalog No	YP-mAb-18822
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CRIP2 CRP2
Immunogen	Synthesized peptide derived from human CRIP2
Specificity	This antibody detects endogenous levels of CRIP2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CRIP;CRP2;ESP1;CRIP2;CRIP2 CRP2;Cysteine-rich protein 2 (CRP-2) (Protein ESP1)
Observed Band	22kDa
Calculated Molecular Weight	22kDa
Cell Pathway	
Tissue Specificity	Widespread tissue expression; highest levels in the heart.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images