



MOCS3 Mouse mAb

Catalog No	YP-mAb-18841
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MOCS3 UBA4
Immunogen	Synthesized peptide derived from human MOCS3
Specificity	This antibody detects endogenous levels of MOCS3 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UBA4;MOCS3;MOCS3 UBA4;Adenylyltransferase and sulfurtransferase MOCS3 (Molybdenum cofactor synthesis protein 3) (Molybdopterin synthase sulfurylase) (MPT synthase sulfurylase) [Includes: Molybdopterin-synthase adenylyltransferase (Adenylyltransferase MOCS3) (Sulfur carrier protein MOCS2A adenylyltransferase);Molybdopterin-synthase sulfurtransferase (Sulfur carrier protein MOCS2A sulfurtransferase) (Sulfurtransferase MOCS3)]
Observed Band	50kDa
Calculated Molecular Weight	50kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	Plays a central role in 2-thiolation of mcm(5)S(2)U at tRNA wobble positions of cytosolic tRNA(Lys), tRNA(Glu) and tRNA(Gln). Also essential during biosynthesis of the molybdenum cofactor. Acts by mediating the C-terminal thiocarboxylation of sulfur carriers URM1 and MOCS2A. Its N-terminus first activates URM1 and MOCS2A as acyl-adenylates (-COAMP), then the persulfide sulfur on the catalytic cysteine is transferred to URM1 and MOCS2A to form thiocarboxylation (-COSH) of their C-terminus. The reaction probably involves hydrogen sulfide that is generated from the persulfide intermediate and



that acts as nucleophile towards URM1 and MOCS2A. Subsequently, a transient disulfide bond is formed. Does not use thiosulfate as sulfur donor; NFS1 probably acting as a sulfur donor for thiocarboxylation reactions.

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images