




CECR5 Mouse mAb

Catalog No	YP-mAb-18852
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 324-423 of human CECR5
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CECR5
Observed Band	46kDa
Calculated Molecular Weight	46kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Predicted to be involved in glycerophospholipid biosynthetic process. Predicted to be active in mitochondrion.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of various lysates using CECR5 Mouse mAb (A20974) at 1:500 dilution. Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 1s.