




HSD17B13 Mouse mAb

Catalog No	YP-mAb-18865
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-211 of human HSD17B13
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLDP;SCDR9;NIIL497;SDR16C3;HMFN0376;HSD17B13
Observed Band	34kDa
Calculated Molecular Weight	34kDa
Cell Pathway	
Tissue Specificity	Liver, kidney, heart, and skeletal muscles. Expressed at moderate levels in the adipose tissue.
Function	
Background	Western blot analysis of extracts of various cell lines, using HSD17B13 antibody (A21056) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25μg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 1s.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of various lysates using KL Mouse pAb (A12028) at 1:3000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25μg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 5s.