



SAMD9L Mouse mAb

Catalog No	YP-mAb-18882
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SAMD9L C7orf6 DRIF2 KIAA2005 UEF
Immunogen	Synthesized peptide derived from human SAMD9L
Specificity	This antibody detects endogenous levels of SAMD9L at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C7orf6;DRIF2;KIAA2005;SAM domain-containing protein 9-like;Sterile alpha motif domain-containing protein 9-like;SAMD9L C7orf6 DRIF2 KIAA2005 UEF; Sterile alpha motif domain-containing protein 9-like (SAM domain-containing protein 9-like)
Observed Band	160-180 kDa
Calculated Molecular Weight	1584 aa, 185 kDa
Cell Pathway	Early endosome .
Tissue Specificity	Widely expressed in adult and fetal tissues. Variable expression in tumors. Down-regulated in breast cancer.
Function	May be involved in endosome fusion. Mediates down-regulation of growth factor signaling via internalization of growth factor receptors.
Background	
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images