




SLIT1 Mouse mAb

Catalog No	YP-mAb-18902
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MEGF4;SLIL1;SLIT3;SLIT-1;SLIT1
Observed Band	168kDa
Calculated Molecular Weight	168kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Enables Roundabout binding activity. Involved in axon extension involved in axon guidance; motor neuron axon guidance; and negative chemotaxis. Predicted to be located in extracellular region. Predicted to be active in extracellular space.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells, using SLIT1 Mouse pAb (A16430) at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3min.