



GPX8 Mouse mAb

Catalog No	YP-mAb-19001
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GPX8 UNQ847/PRO1785
Immunogen	Synthesized peptide derived from human GPX8 AA range: 158-172
Specificity	This antibody detects endogenous levels of GPX8 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GPx-8;UNQ847;EPLA847;GSHPx-8;GPX8;GPX8 UNQ847/PRO1785;Probable glutathione peroxidase 8 (GPx-8) (GSHPx-8) (EC 1.11.1.9)
Observed Band	24kDa
Calculated Molecular Weight	24kDa
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images