



PDIA6 Mouse mAb

Catalog No	YP-mAb-19002
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PDIA6 ERP5 P5 TXNDC7
Immunogen	Synthesized peptide derived from human PDIA6
Specificity	This antibody detects endogenous levels of PDIA6 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	P5;ERP5;TXNDC7;PDIA6;PDIA6 ERP5 P5 TXNDC7;Protein disulfide-isomerase A6 (Endoplasmic reticulum protein 5) (ER protein 5) (ERp5) (Protein disulfide isomerase P5) (Thioredoxin domain-containing protein 7)
Observed Band	48kDa
Calculated Molecular Weight	48kDa
Cell Pathway	Endoplasmic reticulum lumen . Cell membrane . Melanosome . Identified by mass spectrometry in melanosome fractions from stage I to stage IV (PubMed:12643545). .
Tissue Specificity	Expressed in platelets (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images