



MSRB2 Mouse mAb

Catalog No	YP-mAb-19051
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MSRB2 CBS-1 MSRB CGI-131
Immunogen	Synthesized peptide derived from human MSRB2
Specificity	This antibody detects endogenous levels of MSRB2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CBS1;MSRB;PILB;CBS-1;CGI-131;MSRB2;MSRB2 CBS-1 MSRB CGI-131; Methionine-R-sulfoxide reductase B2;mitochondrial (MsrB2) (EC 1.8.4.-); mitochondrial (MsrB2) (EC 1.8.4.-)
Observed Band	20kDa
Calculated Molecular Weight	20kDa
Cell Pathway	Mitochondrion .
Tissue Specificity	Ubiquitous. Detected in retina, ocular ciliary body, skeletal muscle, heart, colon, bone marrow, cerebellum, small intestine, fetal brain, fetal liver, kidney, spinal cord, lung, placenta and prostate.
Function	Methionine-sulfoxide reductase that specifically reduces methionine (R)-sulfoxide back to methionine. While in many cases, methionine oxidation is the result of random oxidation following oxidative stress, methionine oxidation is also a post-translational modification that takes place on specific residue. Upon oxidative stress, may play a role in the preservation of mitochondrial integrity by decreasing the intracellular reactive oxygen species build-up through its scavenging role, hence contributing to cell survival and protein maintenance.
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images