



EXPH5 Mouse mAb

Catalog No	YP-mAb-19139
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EXPH5 KIAA0624 SLAC2B
Immunogen	Synthesized peptide derived from human EXPH5
Specificity	This antibody detects endogenous levels of EXPH5 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DKFZp586F1223;DKFZp781H795;Exophilin 5;Exophilin5;KIAA624;MGC133291; MGC134967;SLAC2-B;SLAC2B;slp homolog lacking C2 domains b; synaptotagmin-like homologue lacking C2 domains b;synaptotagmin-like protein homolog lacking C2 domains b;EXPH5 KIAA0624 SLAC2B;Exophilin-5 (Synaptotagmin-like protein homolog lacking C2 domains b) (Slac2-b) (Slp homolog lacking C2 domains b)
Observed Band	
Calculated Molecular Weight	23 kDa; 223kD(Calculated).
Cell Pathway	
Tissue Specificity	Expressed in keratinocytes.
Function	May act as Rab effector protein and play a role in vesicle trafficking.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images