



CXorf22 Mouse mAb

Catalog No	YP-mAb-19151
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CXorf22
Immunogen	Synthesized peptide derived from human CXorf22
Specificity	This antibody detects endogenous levels of CXorf22 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CXorf22;Uncharacterized protein CXorf22
Observed Band	
Calculated Molecular Weight	351kD
Cell Pathway	Cytoplasm, cytoskeleton, flagellum basal body .
Tissue Specificity	Highly expressed in spermatozoa (at protein level).
Function	Plays a role in flagellar formation and sperm motility.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images