



THADA Mouse mAb

Catalog No	YP-mAb-19154
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	THADA GITA KIAA1767
Immunogen	Synthesized peptide derived from human THADA
Specificity	This antibody detects endogenous levels of THADA at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GITA;ARMC13;Trm732;THADA;THADA GITA KIAA1767;Thyroid adenoma-associated protein (Gene inducing thyroid adenomas protein)
Observed Band	220kDa
Calculated Molecular Weight	220kDa
Cell Pathway	
Tissue Specificity	Expressed in pancreas, adrenal medulla, thyroid, adrenal cortex, testis, thymus, small intestine and stomach.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

