



NA14 Mouse mAb

Catalog No	YP-mAb-19157
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SSNA1 NA14
Immunogen	Synthesized peptide derived from human NA14
Specificity	This antibody detects endogenous levels of NA14 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	N14;NA 14;NA14;Nuclear autoantigen of 14 kDa;Sjogren syndrome nuclear autoantigen 1;Sjogren syndrome nuclear autoantigen 1;Sjogren's syndrome nuclear autoantigen 1;SSNA 1;SSNA1 NA14;Sjogren syndrome nuclear autoantigen 1 (Nuclear autoantigen of 14 kDa)
Observed Band	
Calculated Molecular Weight	14 kDa; 14kD(Calculated).
Cell Pathway	Nucleus. Cytoplasm, cytoskeleton, microtubule organizing center, centrosome .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

