



CHMP1B Mouse mAb

Catalog No	YP-mAb-19161
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CHMP1B C18orf2
Immunogen	Synthesized peptide derived from human CHMP1B
Specificity	This antibody detects endogenous levels of CHMP1B at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Vps46B;C10orf2;C18orf2;CHMP1.5;Vps46-2;C18-ORF2;hVps46-2;CHMP1B;CHMP1B C18orf2;Charged multivesicular body protein 1b (CHMP1.5) (Chromatin-modifying protein 1b) (CHMP1b) (Vacuolar protein sorting-associated protein 46-2) (Vps46-2) (hVps46-2)
Observed Band	30kDa
Calculated Molecular Weight	22kDa
Cell Pathway	Cytoplasm, cytosol. Endosome. Late endosome membrane ; Peripheral membrane protein . Localizes to the midbody of dividing cells, colocalizing with CEP55 and CHMP5. Localized at the periphery of the Fleming body.
Tissue Specificity	Widely expressed. Expressed in pancreas, kidney, skeletal muscle, liver, lung, placenta and brain.
Function	Probable peripherally associated component of the endosomal sorting required for transport complex III (ESCRT-III) which is involved in multivesicular bodies (MVBs) formation and sorting of endosomal cargo proteins into MVBs. MVBs contain intraluminal vesicles (ILVs) that are generated by invagination and scission from the limiting membrane of the endosome and mostly are delivered to lysosomes enabling degradation of membrane proteins, such as stimulated growth factor receptors, lysosomal enzymes and lipids. The MVB pathway



appears to require the sequential function of ESCRT-O, -I, -II and -III complexes. ESCRT-III proteins mostly dissociate from the invaginating membrane before the ILV is released. The ESCRT machinery also functions in topologically equivalent membrane fission events, such as the terminal stages of cytokinesis and the budding of enveloped viruses (HIV-1 and other lentiviruses). ESCRT-III proteins are believed to mediate the necessary vesicle extrusion and/or membrane fission activities, possibly in conjunction with the AAA ATPase VPS4. Involved in cytokinesis. Involved in recruiting VPS4A and/or VPS4B and SPAST to the midbody of dividing cells. Involved in HIV-1 p6- and p9-dependent virus release.

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

