




SRR Mouse mAb

Catalog No	YP-mAb-19173
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ILV1;ISO1;SRR
Observed Band	37kDa
Calculated Molecular Weight	37kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Enables several functions, including L-serine ammonia-lyase activity; PDZ domain binding activity; and anion binding activity. Involved in pyruvate biosynthetic process; response to lipopolysaccharide; and serine family amino acid metabolic process. Located in cytoplasm and neuronal cell body
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of various lysates using SRR Mouse mAb (A9325) at 1:3000 dilution. Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25μg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 60s.