



DDX56 Mouse mAb

Catalog No	YP-mAb-19180
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DDX56 DDX21 NOH61
Immunogen	Synthesized peptide derived from human DDX56
Specificity	This antibody detects endogenous levels of DDX56 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DDX21;DDX26;NOH61;DDX56;DDX56 DDX21 NOH61;Probable ATP-dependent RNA helicase DDX56 (ATP-dependent 61 kDa nucleolar RNA helicase) (DEAD box protein 21) (DEAD box protein 56)
Observed Band	61kDa
Calculated Molecular Weight	62kDa
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Detected in heart, brain, liver, pancreas, placenta and lung.
Function	May play a role in later stages of the processing of the pre-ribosomal particles leading to mature 60S ribosomal subunits. Has intrinsic ATPase activity.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

