



STAMBPL1 Mouse mAb

Catalog No	YP-mAb-19215
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	STAMBPL1 AMSHLP KIAA1373
Immunogen	Synthesized peptide derived from human STAMBPL1
Specificity	This antibody detects endogenous levels of STAMBPL1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AMSH-FP;AMSH-LP;ALMalpha;bA399O19.2;STAMBPL1;STAMBPL1 AMSHLP KIAA1373;AMSH-like protease (AMSH-LP) (STAM-binding protein-like 1)
Observed Band	50kDa
Calculated Molecular Weight	50kDa
Cell Pathway	Ubiquitously expressed.
Tissue Specificity	
Function	Zinc metalloprotease that specifically cleaves 'Lys-63'-linked polyubiquitin chains. Does not cleave 'Lys-48'-linked polyubiquitin chains.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

