



SLC35F2 Mouse mAb

Catalog No	YP-mAb-19234
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC35F2
Immunogen	Synthesized peptide derived from human SLC35F2
Specificity	This antibody detects endogenous levels of SLC35F2 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HSNOV1;SLC35F2;Solute carrier family 35 member F2
Observed Band	41-50 kDa
Calculated Molecular Weight	374 aa, 41 kDa
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	Putative solute transporter.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images