



DNAJC1 Mouse mAb

Catalog No	YP-mAb-19275
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DNAJC1 HTJ1
Immunogen	Synthesized peptide derived from human DNAJC1
Specificity	This antibody detects endogenous levels of DNAJC1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DnaJ protein homolog MTJ1;DNAJC1;DNAJL1;ERdj1;HTJ1;MTJ1;DNAJC1 HTJ1;DnaJ homolog subfamily C member 1 (DnaJ protein homolog MTJ1)
Observed Band	64 kDa
Calculated Molecular Weight	554 aa, 64 kDa
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type I membrane protein . Nucleus membrane; Single-pass type I membrane protein. Microsome membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	May modulate protein synthesis.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images