




METTL9 Mouse mAb

Catalog No	YP-mAb-19317
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FBP interacting repressor;FIR;PUF60;Ro binding protein 1;RoBP1;RoBPI;Siah binding protein 1;Siah BP1;SIAHBP1;PUF60 FIR ROBPI SIAHBP1;Poly(U)-binding-splicing factor PUF60 (60 kDa poly(U)-binding-splicing factor) (FUSE-binding protein-interacting repressor) (FBP-interacting repressor) (Ro-binding protein 1) (RoBP1) (Siah-binding protein 1) (Siah-BP1);DREV;PAP1;DREV1;CGI-81;hMETTL9;METTL9
Observed Band	60 kDa
Calculated Molecular Weight	60 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Enables protein-L-histidine N-pros-methyltransferase activity. Predicted to be involved in methylation. Is active in endoplasmic reticulum.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of various lysates using METTL9 Mouse mAb (A20313) at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25μg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 30s.