



COX8A Mouse mAb

Catalog No	YP-mAb-19329
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TRIM9 KIAA0282 RNF91;E3 ubiquitin-protein ligase TRIM9 (RING finger protein 91) (Tripartite motif-containing protein 9);COX;COX8;VIII;COX8L;COX8-2;VIII-L;MC4DN15;COX8A
Observed Band	8kDa
Calculated Molecular Weight	78kD
Cell Pathway	
Tissue Specificity	
Function	
Background	The protein encoded by this gene is the terminal enzyme of the respiratory chain, coupling the transfer of electrons from cytochrome c to molecular oxygen, with the concomitant production of a proton electrochemical gradient across the inner mitochondrial membrane. In addition to 3 mitochondrially encoded subunits, which perform the catalytic function, the eukaryotic enzyme contains nuclear-encoded smaller subunits, ranging in number from 4 in some organisms to 10 in mammals. It has been proposed that nuclear-encoded subunits may be involved in the modulation of the catalytic function. This gene encodes one of the nuclear-



encoded subunits.


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from mouse skeletal muscle, using COX8A mouse mAb (A3290).