




PTGR3/ZADH2 Mouse mAb

Catalog No	YP-mAb-19346
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TMEM182 UNQ6974/PRO21957;Transmembrane protein 182;PRG-3;ZADH2; PTGR-3;PTGR3/ZADH2
Observed Band	40kDa
Calculated Molecular Weight	25kD
Cell Pathway	
Tissue Specificity	
Function	
Background	Predicted to enable 13-prostaglandin reductase activity. Predicted to be involved in negative regulation of fat cell differentiation. Predicted to be located in peroxisome
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of various lysates using PTGR3/ZADH2 Mouse mAb (A25037) at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution. Lysates / proteins: 25 μ g per lane. Blocking buffer: 3 % nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 20s