




FAM129B Mouse mAb

Catalog No	YP-mAb-19378
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	OC58;MEG-3;C9orf88;FAM129B;MINERVA;bA356B19.6
Observed Band	95kDa
Calculated Molecular Weight	84kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Enables transcription coactivator activity. Involved in several processes, including gonadotropin secretion; positive regulation of transcription regulatory region DNA binding activity; and regulation of cellular macromolecule biosynthetic process. Located in several cellular components, including adherens junction; cytosol; and nucleoplasm.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of various lysates using FAM129B Mouse mAb (A20173) at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 5s.