



GAAP Monoclonal Antibody

Catalog No	YP-mAb-00401
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMBIM4
Immunogen	The antiserum was produced against synthesized peptide derived from human TMBIM4. AA range:121-170
Specificity	GAAP Monoclonal Antibody detects endogenous levels of GAAP protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cirhin interaction protein;CIRIP;CRYBP1;GAAP;HIV EP1;HIVEP1;MBP 1;PRDII BF1;Schnurri 1;ZAS1;TMBIM4;Protein lifeguard 4;LFG4;CGI-119;Golgi anti-apoptotic protein;Protein S1R;Transmembrane BAX inhibitor motif-containing protein 4;Z-protein
Observed Band	23kD
Calculated Molecular Weight	23kD
Cell Pathway	Golgi apparatus membrane ; Multi-pass membrane protein .
Tissue Specificity	Adrenal gland,Blood,Placenta,
Function	similarity:Belongs to the BI1 family.,
Background	similarity:Belongs to the BI1 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

