



Mammaglobin A Monoclonal Antibody

Catalog No	YP-mAb-00439
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SCGB2A2
Immunogen	The antiserum was produced against synthesized peptide derived from human Mammaglobin. AA range:41-90
Specificity	Mammaglobin A Monoclonal Antibody detects endogenous levels of Mammaglobin A protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	hMAM;Mammaglobin 1;Mammaglobin A;MGB1;SCGB2A2;UGB2;Mammaglobin; Mammaglobin-A;Mammaglobin-1;Secretoglobin family 2A member 2
Observed Band	20kD
Calculated Molecular Weight	20kD
Cell Pathway	Secreted .
Tissue Specificity	Mammary gland specific. Over-expressed in breast cancer.
Function	similarity:Belongs to the secretoglobin family. Lipophilin subfamily.,tissue specificity:Mammary gland specific. Over-expressed in breast cancer.,
Background	similarity:Belongs to the secretoglobin family. Lipophilin subfamily.,tissue specificity:Mammary gland specific. Over-expressed in breast cancer.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

