



MLH3 Monoclonal Antibody

Catalog No	YP-mAb-00451
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MLH3
Immunogen	The antiserum was produced against synthesized peptide derived from human MLH3. AA range:521-570
Specificity	MLH3 Monoclonal Antibody detects endogenous levels of MLH3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HNPCC7;MLH3;mutL homolog 3 (E. coli);MutL protein homolog 3;DNA mismatch repair protein Mlh3
Observed Band	164kD
Calculated Molecular Weight	164kD
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitous.
Function	disease:Defects in MLH3 are a cause of somatic colorectal cancer (CRC) [MIM:114500].,disease:Defects in MLH3 are the cause of hereditary non-polyposis colorectal cancer type 7 (HNPCC7) [MIM:604395]. Mutations in more than one gene locus can be involved alone or in combination in the production of the HNPCC phenotype (also called Lynch syndrome). Most families with clinically recognized HNPCC have mutations in either MLH1 or MSH2 genes. HNPCC is an autosomal, dominantly inherited disease associated with marked increase in cancer susceptibility. It is characterized by a familial predisposition to early onset colorectal carcinoma (CRC) and extra-colonic cancers of the gastrointestinal, urological and female reproductive tracts. HNPCC is reported to be the most common form of inherited colorectal cancer in the Western world,



and accounts for 15% of all colon cancers. Cancers in HNPCC origina

Background

This gene is a member of the MutL-homolog (MLH) family of DNA mismatch repair (MMR) genes. MLH genes are implicated in maintaining genomic integrity during DNA replication and after meiotic recombination. The protein encoded by this gene functions as a heterodimer with other family members. Somatic mutations in this gene frequently occur in tumors exhibiting microsatellite instability, and germline mutations have been linked to hereditary nonpolyposis colorectal cancer type 7 (HNPCC7). Several alternatively spliced transcript variants have been identified, but the full-length nature of only two transcript variants has been determined. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

