



Siah-1/2 Monoclonal Antibody

Catalog No	YP-mAb-00522
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SIAH1/SIAH2
Immunogen	The antiserum was produced against synthesized peptide derived from human SIAH1. AA range:181-230
Specificity	Siah-1/2 Monoclonal Antibody detects endogenous levels of Siah-1/2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	hSIAH1;HUMSIAH;Seven in absentia homolog 1;Siah 1;Siah 1a;SIAH1;Siah-1/2;SIAH1/SIAH2;E3 ubiquitin-protein ligase SIAH1/2;E3 ubiquitin-protein ligase SIAH1;Siah-1;Siah-1a;SIAH2;E3 ubiquitin-protein ligase SIAH2;Seven in absentia homolog 2;Siah-2;hSiah2
Observed Band	34kD
Calculated Molecular Weight	34kD
Cell Pathway	Cytoplasm. Nucleus. Predominantly cytoplasmic. Partially nuclear.
Tissue Specificity	Widely expressed at a low level. Down-regulated in advanced hepatocellular carcinomas.
Function	domain:The RING-type zinc finger domain is essential for ubiquitin ligase activity.,domain:The SBD domain (substrate-binding domain) mediates the homodimerization and the interaction with substrate proteins. It is related to the TRAF family.,function:E3 ubiquitin-protein ligase that mediates ubiquitination and subsequent proteasomal degradation of target proteins. E3 ubiquitin ligases accept ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates. Mediates E3 ubiquitin ligase activity either through direct binding to substrates or



by functioning as the essential RING domain subunit of larger E3 complexes. Triggers the ubiquitin-mediated degradation of many substrates, including proteins involved in transcription regulation (MYB, POU2AF1, PML and RBBP8), a cell surface receptor (DCC), cytoplasmic signal

Background

This gene encodes a protein that is a member of the seven in absentia homolog (SIAH) family. The protein is an E3 ligase and is involved in ubiquitination and proteasome-mediated degradation of specific proteins. The activity of this ubiquitin ligase has been implicated in the development of certain forms of Parkinson's disease, the regulation of the cellular response to hypoxia and induction of apoptosis. Alternative splicing results in several additional transcript variants, some encoding different isoforms and others that have not been fully characterized. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Siah-1/2 Monoclonal Antibody