



UROC28 Monoclonal Antibody

Catalog No	YP-mAb-00547
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PBOV1
Immunogen	The antiserum was produced against synthesized peptide derived from human PBOV1. AA range:31-80
Specificity	UROC28 Monoclonal Antibody detects endogenous levels of UROC28 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UROC28;PBOV1;Prostate and breast cancer overexpressed gene 1 protein; Protein UROC28;UC28
Observed Band	19kD
Calculated Molecular Weight	19kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Expressed in colon, prostate, small intestine, testis and spleen, with lower expression in thymus, ovary, and peripheral blood leukocytes. Up-regulated expression in prostate, breast, and bladder cancer, but not in lung and colon cancer.
Function	tissue specificity:Expressed in colon, prostate, small intestine, testis and spleen, with lower expression in thymus, ovary, and peripheral blood leukocytes. Up-regulated expression in prostate, breast, and bladder cancer, but not in lung and colon cancer.,
Background	This intronless gene encodes a protein of unknown function. Its expression is up-regulated in some types of cancer, including prostate, breast, and bladder cancer. [provided by RefSeq, Aug 2011],



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

