



BCAS4 Monoclonal Antibody

Catalog No	YP-mAb-00682
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BCAS4
Immunogen	The antiserum was produced against synthesized peptide derived from human BCAS4. AA range:31-80
Specificity	BCAS4 Monoclonal Antibody detects endogenous levels of BCAS4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BCAS4;Breast carcinoma-amplified sequence 4
Observed Band	23kD
Calculated Molecular Weight	23kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Brain, thymus, spleen, kidney and placenta. Overexpressed in most breast cancer cell lines.
Function	disease:A chromosomal aberration involving BCAS4 may be a cause of breast cancer. Translocation t(17;20)(q23;q13) with BCAS3.,similarity:Belongs to the cappuccino family.,tissue specificity:Brain, thymus, spleen, kidney and placenta. Overexpressed in most breast cancer cell lines.,
Background	disease:A chromosomal aberration involving BCAS4 may be a cause of breast cancer. Translocation t(17;20)(q23;q13) with BCAS3.,similarity:Belongs to the cappuccino family.,tissue specificity:Brain, thymus, spleen, kidney and placenta. Overexpressed in most breast cancer cell lines.,



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

