



p115 Monoclonal Antibody

Catalog No	YP-mAb-00713
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	USO1
Immunogen	The antiserum was produced against synthesized peptide derived from human USO1. AA range:913-962
Specificity	p115 Monoclonal Antibody detects endogenous levels of p115 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RhoGAP4;Rho-GAP hematopoietic protein C1;Rho GTPase-activating protein 4; RGC1;p115;USO1;General vesicular transport factor p115;VDP;Protein USO1 homolog;Transcytosis-associated protein;TAP;Vesicle-docking protein
Observed Band	108kD
Calculated Molecular Weight	108kD
Cell Pathway	Cytoplasm, cytosol . Golgi apparatus membrane ; Peripheral membrane protein . Recycles between the cytosol and the Golgi apparatus during interphase. During interphase, the phosphorylated form is found exclusively in cytosol; the unphosphorylated form is associated with Golgi apparatus membranes. .
Tissue Specificity	Epithelium,Human skeletal muscle,Liver,Uterus,
Function	domain:Composed of a globular head, an elongated tail (coiled-coil) and a highly acidic C-terminal domain.,function:General vesicular transport factor required for intercisternal transport in the Golgi stack; it is required for transcytotic fusion and/or subsequent binding of the vesicles to the target membrane. May well act as a vesicular anchor by interacting with the target membrane and holding the vesicular and target membranes in proximity.,PTM:Phosphorylated in a cell cycle-specific manner; phosphorylated in interphase but not in mitotic cells.



Dephosphorylated protein associates with the Golgi membrane; phosphorylation promotes dissociation.,similarity:Belongs to the VDP/USO1/EDE1 family., subcellular location:Recycles between the cytosol and the Golgi apparatus during interphase.,

Background

The protein encoded by this gene is a peripheral membrane protein which recycles between the cytosol and the Golgi apparatus during interphase. It is regulated by phosphorylation: dephosphorylated protein associates with the Golgi membrane and dissociates from the membrane upon phosphorylation. Ras-associated protein 1 recruits this protein to coat protein complex II (COPII) vesicles during budding from the endoplasmic reticulum, where it interacts with a set of COPII vesicle-associated SNAREs to form a cis-SNARE complex that promotes targeting to the Golgi apparatus. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2014],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using p115 Monoclonal Antibody