



SFT2D2 Monoclonal Antibody

Catalog No	YP-mAb-00723
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SFT2D2
Immunogen	The antiserum was produced against synthesized peptide derived from human SFT2B. AA range:61-110
Specificity	SFT2D2 Monoclonal Antibody detects endogenous levels of SFT2D2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SFT2 domain-containing protein 2;UNQ512/PRO1027;Vesicle transport protein SFT2B;SFT2D2
Observed Band	18kD
Calculated Molecular Weight	18kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Lymph,
Function	function:May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex.,similarity:Belongs to the SFT2 family.,
Background	function:May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex.,similarity:Belongs to the SFT2 family.,
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

