



ASCL1 Monoclonal Antibody

Catalog No	YP-mAb-01541
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ASCL1
Immunogen	The antiserum was produced against synthesized peptide derived from human ASCL1. AA range:109-158
Specificity	ASCL1 Monoclonal Antibody detects endogenous levels of ASCL1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Achaete-scute homolog 1;ASH1;ASH-1;bHLHa46;Class A basic helix-loop-helix protein 46;ASCL1;HASH1
Observed Band	25kD
Calculated Molecular Weight	25kD
Cell Pathway	Nucleus .
Tissue Specificity	Lung,Thyroid carcinoma,
Function	function:May play a role at early stages of development of specific neural lineages in most regions of the CNS, and of several lineages in the PNS. Essential for the generation of olfactory and autonomic neurons. Activates transcription by binding to the E box (5'-CANNTG-3').,similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein. Forms a heterodimer with TCF3.,
Background	achaete-scute family bHLH transcription factor 1(ASCL1) Homo sapiens This gene encodes a member of the basic helix-loop-helix (BHLH) family of transcription factors. The protein activates transcription by binding to the E box (5'-CANNTG-3'). Dimerization with other BHLH proteins is required



for efficient DNA binding. This protein plays a role in the neuronal commitment and differentiation and in the generation of olfactory and autonomic neurons. Mutations in this gene may contribute to the congenital central hypoventilation syndrome (CCHS) phenotype in rare cases. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ASCL1 Monoclonal Antibody